





Synchrotron Data Workflow AstraZeneca SE



Patrik Johansson



Historical overview 2004-2023




- 3717 structures determined from in-house x-ray or synchrotron
- 22972 X-ray datasets collected at synchrotrons, including 20079 which did not result in a structure




Figures adapted from *J Synchrotron Radiat.* 2025 Feb 6;32(Pt 2):294–303



Synchrotrons play a critical role in our workflow



- Near-weekly synchrotron trips




- Using multiple synchrotrons to ensure coverage around the year


Figures adapted from *J Synchrotron Radiat.* 2025 Feb 6;32(Pt 2):294–303



Experiment setup & metadata export



Data import and processing



Finding the files

LUND

`$data_dir/$tripid*/process*/$alias/$alias-$crysid/xds_*/MAXIV*/Edna*/results/*final.mtz`

`$data_dir/$tripid*/process*/$alias/$alias-$crysid/xds_*/MAXIV*/AutoPROCTask_*/results/*staraniso_alldata-unique.mtz`

..

ESRF

`$data_dir/$tripid*/PROCESSED_DATA*/$alias/$alias*-$crysid/run*/autoprocessing*/EDNA*/results/*_anom_truncate.mtz`

`$data_dir/$tripid*/PROCESSED_DATA*/$alias/$alias*$crysid/run*/**/autoPROC/autoPROC_staraniso/*staraniso_alldata-unique.mtz`

..

Diamond

`$data_dir/$tripid*/processed/$alias/*$crysid/*$crysid*/**/autoPROC/staraniso_alldata-unique.mtz`

`$data_dir/$tripid*/processed/$alias/*$crysid/*$crysid*/**/xia2*/DataFiles/*free.mtz"`

..

SLS

`$data_dir/$tripid*/$alias/$alias-$crysid*/adp_*/autoproc/*/*staraniso_alldata-unique.mtz`

..


Soleil

`$data_dir/$tripid*/PROCESSED_DATA/$alias/$alias-$crysid/xdsme_*/ccp4if/$alias-$crysid*.mtz`

..

+ associated logfiles

Data import and processing



Summary, synchrotron data/logistics

- Common ISPYB / EXI architecture very helpful for end users
 - Variations in the details of input / output a complicating factor
- Data retrieval different for every synchrotron
 - Requirements from corporate IT sometimes not compatible with available options
 - DART download, managed by in-house IT
- Remote access
 - Direct access not compatible with IT safety rules
 - MXCUBE3
- Changes of location and names of preprocessed files / logs
 - Requires adaptation of in-house pipeline

